

■ Reagenzglasgestell mit Farbkennung

Kennzeichnung durch Farbköpfe (weiß, orange, gelb, blau, grün, rot)
1 Satz Farbköpfe inklusive
Autoklavierbar

■ Test tube rack with colour identification


Identification by coloured knobs (white, orange, yellow, blue, green, red)
1 set of coloured knobs included
Autoclavable

■ Porte-tubes à essai avec identification de couleur

Identification par boutons colorés (blanc, orange, jaune, bleue, vert, rouge)
1 set de boutons inclus
Autoclavable

■ Gradillas para tubos de ensayos con identificación de color

Identificación por botones de color (blanco, naranja, amarillo, azul, verde, rojo)
1 kit de botones incluido
Autoclavable




Material	A	B	C	D	Me/Qty	Me/Qty	Artikel Nr.
18/10-STAHL	248	104	60	13	5x12	60	8250
18/10-STAHL	248	104	85	16	5x12	60	8251
18/10-STAHL	248	104	85	20	4x10	40	8252
18/10-STAHL	248	104	85	25	3x7	21	8253
18/10-STAHL	248	104	100	30	3x7	21	8254

■ Reagenzglasgestell

■ Test tube stand

■ Porte-tubes à essai

■ Gradillas para tubos de ensayo




Material	A	B	C	D	Typ	Artikel Nr.
18/10-STAHL	145	52	90	17	2x6	8210
18/10-STAHL	278	52	90	17	2x12	8212
18/10-STAHL	287	97	90	17	4x12	8214
18/10-STAHL	164	60	90	20	2x6	8220
18/10-STAHL	310	60	90	20	2x12	8222

■ Reagenzglasgestell, rund

■ Test tube stand, round

■ Porte tubes à essai, forme ronde

■ Gradillas redondas para tubos de ensayo



Material	A	C	D	Typ	Artikel Nr.
18/10-STAHL	98	165	20	9 x	8230
18/10-STAHL	116	165	20	12x	8231
18/10-STAHL	194	165	20	45x	8233

■ Reagenzglasgestell Z-Form

■ Test tube stand Z-shape

■ Porte-tubes à essai, forme Z

■ Gradillas para tubos de ensayo, forma Z



Material	A	B	C	D	Typ	Artikel Nr.
18/10-STAHL	50	165	85	13	10x5	8290
18/10-STAHL	50	220	110	17	10x5	8291
18/10-STAHL	50	85	85	13	5x5	8292
18/10-STAHL	50	110	110	17	5x5	8293

■ **Transportwagen für Transportkörbe 18/10 Stahl**

Mit Griff, 4 Lenkrollen (davon 2 feststellbar)

■ **Transport cart for transport baskets 18/10 stainless steel**


With handle, 4 guide rolls (2 lockable)

■ **Chariot pour paniers de transport acier inox 18/10**

Avec manche, 4 roues de guidage (2 avec treins)

■ **Carro para cestillos de transporte 18/10 acero inox**

Con mango, 4 ruedas (2 con freno)



Material	A	B	C	Artikel Nr.
18/10 E-POLI	603	403	950	17330

■ **Deckel für Transportkorb 18/10 Stahl**

3 verschiedene Maschenweiten (m=10x10x3, 20x20x3, 35x35x3 mm)

■ **Cover for transport basket 18/10 stainless steel**


3 different mesh dimensions (m=10x10x3, 20x20x3, 35x35x3 mm)

■ **Couvercles pour paniers de transport acier inox 18/10**

3 différents dimensions de mailles (m=10x10x3, 20x20x3, 35x35x3 mm)

■ **Tapas para cestillos de transporte 18/10 acero inox**

3 diferentes tamaños de malla (m=10x10x3, 20x20x3, 35x35x3 mm)



Material	A	B	m	Artikel Nr.
18/10 E-POLI	585	385	10x10x3	17315
18/10 E-POLI	285	385	10x10x3	17316
18/10 E-POLI	585	385	20x20x3	17320
18/10 E-POLI	285	385	20x20x3	17321
18/10 E-POLI	585	385	35x35x3	17325
18/10 E-POLI	285	385	35x35x3	17326

■ **Abtropfgestell mit Abtropfschale**

Für Reagenzgläser und Kolben

■ **Draining rack with drip pan**


For test tubes and flasks

■ **Egouttoir avec cuvette**

Pour tubes à essai et ballons

■ **Panel de escurrido de sobremesa**

Con bandeja de escurrido. Para material de vidrio diverso




Material	A	B	C	Me/Qty	Artikel Nr.
PA-COATING	420	160	610	65	10130
PA-COATING	420	160	610	44	10140

■ **Reagenzglasgestelle**

■ **Test tube stands**

■ **Portes-tubes à essai**

■ **Gradilla para tubos de ensayo**



Material	A	B	C	Me/Qty	Me/Qty D	Artikel Nr.
18/10 E-POLI	130	50	70	2x6	12 16	10170
18/10 E-POLI	210	50	70	2x10	20 16	10171
18/10 E-POLI	250	50	70	2x12	24 16	10172
18/10 E-POLI	250	70	70	3x12	36 16	10173
18/10 E-POLI	130	90	70	4x6	24 16	10174
18/10 E-POLI	250	90	70	4x12	48 16	10175
18/10 E-POLI	210	210	70	10x10	100 16	10176